

## LABORATORY INSTRUCTIONS

### Vacutainer Tube - Order of Draw

Saint Francis Hospital  
 6161 South Yale Avenue  
 Tulsa, OK 74136  
 (918) 494-1300  
 (918) 494-1399 fax

COLOR	COLLECTION TUBE	MIX BY INVERTING
	Blood Cultures	
	Citrate	8 to 10 times
	Clot Activator - 6 ml	5 times
	SST Gel Separator	5 times
	Heparin	8 to 10 times
	EDTA	8 to 10 times
	EDTA	8 to 10 times